

# CYBERKNIFE® STARTUP KIT\*

For your convenience Accuray is pleased to offer you the Startup Kit for the CyberKnife® System. The following ancillary equipment is required to commission or sustain proper functionality of your system. The Startup Kit can be user configured by selecting one or all of the items below. You may also choose to purchase those items on your own\*\*. Please contact Accuray Physics Support with technical questions and your Accuray Sales Person for procurement.

**Physics and QA Equipment** to commission and periodically validate your CyberKnife and treatment planning systems.




Choice of **Water Tanks** to commission the CyberKnife System.



**SF6 Gas Cylinder** for linear accelerator.



**Xsight® Lung Tracking Phantom kit** to QA the Xsight Lung Tracking System.



**In Room Video Communication Set** to continuously communicate with the patient from the control console.



**Laser Kit** to locate and align patient prior to treatment.



Item	Description	Part Number	Vendor
<b>Physics and QA Equipment</b>	SuperMAX Electrometer, 20 meter, Extension Cable [60'] with Tri-axial BNC or TNC, Ion Chambers: A16 Exradin Ion Chamber and A12S Exradin Waterproof Modified Farmer Type Chamber, with Cobalt Calibration Point in Water, Stereotactic Dose Verification Phantom with SRS Dosimetric QA Slab and Plugs for A16 and A12S Ion Chambers, Survey Meter, Pressurized Red/YL, Epson Expression 10000XL Color Flatbed Photo Scanner 2400 dpi Resolution, Electron/Mass Density Phantom - 19 Plugs, GafChromic EBT2 Film, 8" x 10" Sheets, 25 Sheets/Box, Precision Aneroid Barometer 870-1050 hPa w/ Calibration, Glass Mercury Thermometer with 0.2 C Divisions, QCKV-1 Phantom with CyberKnife Positioning Stand, Unfors Solo RF	1010006	Standard Imaging
<b>Choice of Water Tanks</b>	M3000RS Water Tank + 2 Diodes	1013391	ARM
	MP3-M Water Phantom, TPR Option for MP3-M/MP3-S, (2) Dosimetry Diodes	1012236	PTW
<b>SF6 Gas Cylinder</b>	SF6, Sulfur Hexafluoride 99.9% Q, SF6, **Available in the U.S. ONLY**	1010008	Praxair
<b>Xsight® Lung Tracking Phantom kit</b>	Phantom Torso Xsight Lung Tracking, Phantom Cylinder Xsight Lung Tracking, Phantom Ball Cube Xsight Lung Tracking, Phantom Baseplate, Phantom LED Platform, ASSY Xsight Lung Tracking phantom	025190	CIRS
<b>In Room Video Communication Set</b>	Everfocus Keyboard Controller, (4) Everfocus PTZ Cameras, (4) Everfocus Ceiling Mount Adapters, 4-C Color Quad Processor, Everfocus 17" LCD Monitor, AC Power Supply, Siamese Video/Power Cable, 18/2 cm[1K SPOOL], 15 RG-59 BNC Twist on Connectors, Aiphone Master Station, Aiphone Door Station, Aiphone Power Supply **Installation included**	1010007	Cytech
<b>Laser Kit</b>	LAP Red Astor Sagittal Line Laser, (2) LAP Red Astor Crosshair Lasers, Laser Installation on "Zero Vibration" Mount Surfaces **Installation included**	1010005	LAP



Worldwide Headquarters 1310 Chesapeake Terrace | Sunnyvale, CA 94089 | Tel: +1.408.716.4600 | www accuray.com

\*Products listed are not designed, tested or manufactured by Accuray nor does Accuray take any responsibility for supporting these products. All warranties at the warranties of the original manufacturers.  
\*\*Please note that the startup kit should be ordered at least 10 weeks in advance of CyberKnife System acceptance.