

STANDARD *TomoTherapy*[®] QUALITY ASSURANCE KIT

For your convenience Accuray is pleased to offer you the Quality Assurance Kit for the TomoTherapy[®] System. The following ancillary equipment is required to commission or sustain proper functionality of your system. The Quality Assurance Kit can be user configured by selecting one or more of the items below. You may also choose to purchase those items on your own. Please contact Accuray Physics Support with technical questions and your Accuray Sales Person for procurement.

Standard TomoTherapy[®] Quality Assurance Kit is a standard package which includes quality monitoring and performance testing phantoms.

Cheese Phantom



Ion Chambers



Density Plugs



Water Phantoms



Water Tank Kit for obtaining and analyzing radiation beam data.



TomoElectrometer


Film Scanner Kit is used to digitize films for patient-specific quality assurance.



Electrometer for acquiring and analyzing ion chamber data.



Film Analysis Kit for quality assurance and commissioning procedures.



| Item | Description | Part Number |
|--|--|-----------------|
| Standard Quality Assurance Kit | <p>Standard TomoTherapy Quality Assurance package includes: Cheese Phantom, solid water phantom, and ion chamber</p> <p>Additional details:</p> <p>Ion chambers Ion Chamber, 2A1SL, ExradinKit, Exradin A17 Ion Chamber,</p> <p>Virtual Water Phantom 55x15x5cm, no holes 55x15x5cm, drilled cavity at 1cm depth for A1SL Ion Chamber, 55x15x2cm, no holes 55x15x2cm, two drilled cavities at 1cm depth for A1SL ion chamber, 55x15x1cm, no holes 55x15x0.5cm, no holes</p> <p>Cheese Phantom with Density Plugs for Patient QA TomoPhantom (Cheese Phantom) Kit, plugs, contrast, TomoPhantom</p> | 104508-000 |
| Water Tank Kit | <p>TomoScanner Water Tank TomoElectrometer w/ADCL Calibration Computer and TEMS software</p> | H-1006-0000-000 |
| Electrometer Kit <i>**Needed if you did not order Water Tank kit</i> | <p>TomoElectrometer w/ADCL Calibration Computer and TEMS software`</p> | 104940 |
| Film Scanner Kit | <p>Film Digitizer Kit - with Vidar DosimetryPRO Red Scanner Computer Film analysis software</p> | H-1005-0001-000 |
| Film Analysis Kit <i>**If choosing alternative scanner</i> | <p>PC and Software Kit for Existing Vidar DosimetryPRO Advantage or Flatbed Computer Film analysis software</p> | H-1004-0001-000 |



Worldwide Headquarters 1310 Chesapeake Terrace | Sunnyvale, CA 94089 | Tel: +1.408.716.4600 | www accuray.com

*Products listed are not designed, tested or manufactured by Accuray nor does Accuray take any responsibility for supporting these products. All warranties at the warranties of the original manufacturers.
**Please note that the startup kit should be ordered at least 10 weeks in advance of Tomotherapy System acceptance.